

产品名称: DDAH1 Rabbit Monoclonal Antibody
产品货号: AMRe85497



产品概述 (Summary)

产品名称 (Production Name)	DDAH1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DDAH1
别名 (Alternative Names)	DDAH; DDAHI; DDAH-1; HEL-S-16
基因 ID (Gene ID)	23576.0
蛋白 ID (SwissProt ID)	O94760.A synthetic peptide of human DDAH1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 31 kDa; Observed MW: 31 kDa

研究背景 (Background)

Hydrolyzes N(G),N(G)-dimethyl-L-arginine (ADMA) and N(G)-monomethyl-L-arginine (MMA) which act as inhibitors of


产品名称: DDAH1 Rabbit Monoclonal Antibody
产品货号: AMRe85497



NOS. Has therefore a role in the regulation of nitric oxide generation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of DDAH1 in HeLa lysates using DDAH1 antibody.

注意事项 (Note)

For research use only .