

产品名称: DDX5 Rabbit Monoclonal Antibody
产品货号: AMRe85499



产品概述 (Summary)

产品名称 (Production Name)	DDX5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DDX5
别名 (Alternative Names)	DDX5; G17P1; HELR; HLR1; Probable ATP-dependent RNA helicase DDX5; DEAD box protein 5; RNA helicase p68
基因 ID (Gene ID)	1655.0
蛋白 ID (SwissProt ID)	P17844.A synthetic peptide of human DDX5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 69 kDa; Observed MW: 69 kDa

研究背景 (Background)


产品名称: DDX5 Rabbit Monoclonal Antibody
产品货号: AMRe85499




RNA-dependent ATPase activity. The rate of ATP hydrolysis is highly stimulated by single-stranded RNA. May be involved in pre-mRNA splicing.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of DDX5 in C6, 3T3, HeLa lysates using DDX5 antibody.



Immunocytochemistry analysis of DDX5 (green) in HeLa using DDX5 antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using DDX5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: DDX5 Rabbit Monoclonal Antibody
产品货号: AMRe85499



注意事项 (Note)

For research use only .