

产品名称: DYNLL1 Rabbit Monoclonal Antibody
产品货号: AMRe85520



产品概述 (Summary)

产品名称 (Production Name)	DYNLL1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DYNLL1
别名 (Alternative Names)	LC8; PIN; DLC1; DLC8; LC8a; DNCL1; hdlc1; DNCLC1
基因 ID (Gene ID)	8655.0
蛋白 ID (SwissProt ID)	P63167.A synthetic peptide of human DYNLL1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 10 kDa; Observed MW: 10 kDa

研究背景 (Background)

Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be

产品名称: DYNLL1 Rabbit Monoclonal Antibody
产品货号: AMRe85520




involved in linking dynein to cargos and to adapter proteins that regulate dynein function. Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules. May play a role in changing or maintaining the spatial distribution of cytoskeletal structures.


研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Western blot analysis of DYNLL1 in 3T3, HeLa lysates using DYNLL1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using DYNLL1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .