

产品名称: EHD1 Rabbit Monoclonal Antibody
产品货号: AMRe85524



产品概述 (Summary)

产品名称 (Production Name)	EHD1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EHD1
别名 (Alternative Names)	EHD1; H PAST; hPAST1; PAST; PAST homolog 1; PAST1; Tilin
基因 ID (Gene ID)	10938.0
蛋白 ID (SwissProt ID)	Q9H4M9.A synthetic peptide of human EHD1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 61 kDa; Observed MW: 61 kDa

研究背景 (Background)


Acts in early endocytic membrane fusion and membrane trafficking of recycling endosomes.

产品名称: EHD1 Rabbit Monoclonal Antibody
产品货号: AMRe85524




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of EHD1 in K562, rat Brain, C6, 3T3, HeLa lysates using EHD1 antibody.



Immunocytochemistry analysis of EHD1 (green) in HeLa using EHD1 antibody, and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using EHD1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: EHD1 Rabbit Monoclonal Antibody
产品货号: AMRe85524



注意事项 (Note)

For research use only .