

产品名称: eIF1A Rabbit Monoclonal Antibody
产品货号: AMRe85525



产品概述 (Summary)

产品名称 (Production Name)	eIF1A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	eIF1A
别名 (Alternative Names)	EIF1A; EIF4C; eIF-1A; eIF-4C; EIF1AP1
基因 ID (Gene ID)	1964.0
蛋白 ID (SwissProt ID)	P47813.A synthetic peptide of human eIF1A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 16 kDa; Observed MW: 16 kDa

研究背景 (Background)

Seems to be required for maximal rate of protein biosynthesis. Enhances ribosome dissociation into subunits and stabilizes


产品名称: eIF1A Rabbit Monoclonal Antibody
产品货号: AMRe85525




the binding of the initiator Met-tRNA(I) to 40 S ribosomal subunits.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of eIF1A in mouse brain, rat brain lysates using eIF1A antibody.




Western blot analysis of eIF1A in LnCap, K562, U251MG, HepG2, MCF-7, Hela lysates using eIF1A antibody



Immunocytochemistry analysis of eIF1A (green) in 293 using eIF1A antibody, and DAPI (blue).

产品名称: eIF1A Rabbit Monoclonal Antibody
产品货号: AMRe85525



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using eIF1A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .