

产品名称: eIF4A3 Rabbit Monoclonal Antibody
产品货号: AMRe85530



产品概述 (Summary)

产品名称 (Production Name)	eIF4A3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	eIF4A3
别名 (Alternative Names)	DDX48; EIF4A3; eIF4AIII; hNMP 265; NMP 265; NMP265; NUK34
基因 ID (Gene ID)	9775.0
蛋白 ID (SwissProt ID)	P38919.A synthetic peptide of human Eif4a3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)

ATP-dependent RNA helicase. Component of a splicing-dependent multiprotein exon junction complex (EJC) deposited at


产品名称: eIF4A3 Rabbit Monoclonal Antibody
产品货号: AMRe85530




splice junction on mRNAs. The EJC is a dynamic structure consisting of a few core proteins and several more peripheral nuclear and cytoplasmic associated factors that join the complex only transiently either during EJC assembly or during subsequent mRNA metabolism.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of eIF4A3 in HeLa, CHO-K1, C6 lysates using eIF4A3 antibody.



Immunocytochemistry analysis of Eif4a3 (green) in HeLa using Eif4a3 antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Eif4a3 antibody. High-pressure and temperature

产品名称: eIF4A3 Rabbit Monoclonal Antibody
产品货号: AMRe85530



Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .