

产品名称: Emerin Rabbit Monoclonal Antibody
产品货号: AMRe85536



产品概述 (Summary)

产品名称 (Production Name)	Emerin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Emerin
别名 (Alternative Names)	STA; EDMD; LEMD5
基因 ID (Gene ID)	2010.0
蛋白 ID (SwissProt ID)	P50402.A synthetic peptide of human Emerin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 35 kDa

研究背景 (Background)

Emerin may regulate gene expression through binding to other transcriptional regulators. Emerin binds to β -catenin and


产品名称: Emerin Rabbit Monoclonal Antibody
产品货号: AMRe85536




inhibits its nuclear accumulation. Recent studies demonstrate that Emerin is required for HIV-1 infectivity.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Emerin in 3T3, HeLa lysates using Emerin antibody.



Immunocytochemistry analysis of Emerin(green) in HeLa using Emerin antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Emerin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Emerin Rabbit Monoclonal Antibody
产品货号: AMRe85536



注意事项 (Note)

For research use only .