

产品名称: EPS15R Rabbit Monoclonal Antibody
产品货号: AMRe85543



产品概述 (Summary)

产品名称 (Production Name)	EPS15R Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EPS15R
别名 (Alternative Names)	EPS15R
基因 ID (Gene ID)	58513.0
蛋白 ID (SwissProt ID)	Q9UBC2.A synthetic peptide of human EPS15R

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 94 kDa; Observed MW: 130 kDa

研究背景 (Background)

Seems to be a constitutive component of clathrin-coated pits that is required for receptor-mediated endocytosis. Involved


产品名称: EPS15R Rabbit Monoclonal Antibody
产品货号: AMRe85543



in endocytosis of integrin beta-1 (ITGB1) and transferrin receptor (TFR); internalization of ITGB1 as DAB2-dependent cargo but not TFR seems to require association with DAB2.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of EPS15R in K562 lysates using EPS15R antibody.

注意事项 (Note)

For research use only .