

产品名称: Estrogen Receptor alpha Rabbit Monoclonal Antibody
产品货号: AMRe85547



产品概述 (Summary)

产品名称 (Production Name)	Estrogen Receptor alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Estrogen Receptor alpha
别名 (Alternative Names)	ESR1; Era; Eralpha; Estrogen receptor; Estradiol receptor; ER-alpha; Estrogen receptor 1; NR3A1; ER; ESR; ESRA; Estrogen receptor alpha
基因 ID (Gene ID)	2099.0
蛋白 ID (SwissProt ID)	P03372.A synthetic peptide of human Estrogen Receptor alpha

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 66 kDa; Observed MW: 66 kDa

产品名称: Estrogen Receptor alpha Rabbit Monoclonal Antibody
产品货号: AMRe85547




研究背景 (Background)


ER (estrogen receptor 1) a member of the steroid receptor superfamily, contains highly conserved DNA binding (DBD) and ligand binding domains (LBD). Through its estrogen-independent and estrogen-dependent activation domains (AF-1 and AF-2, respectively), ER regulates transcription by recruiting coactivator proteins and interacting with general transcriptional machinery. Phosphorylation provides an important mechanism to regulate ER activity. ER is phosphorylated on multiple sites.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of Estrogen Receptor alpha in K562, HeLa lysates using Estrogen Receptor alpha antibody.



Immunocytochemistry analysis of Estrogen Receptor alpha (green) in MCF-7 using Estrogen Receptor alpha antibody, and DAPI(blue)

产品名称: Estrogen Receptor alpha Rabbit Monoclonal Antibody

产品货号: AMRe85547



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Estrogen Receptor alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .