

产品名称: FADD Rabbit Monoclonal Antibody
产品货号: AMRe85554



产品概述 (Summary)

产品名称 (Production Name)	FADD Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FADD
别名 (Alternative Names)	Mort1/FADD
基因 ID (Gene ID)	14082.0
蛋白 ID (SwissProt ID)	Q61160.A synthetic peptide of mouse FADD

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 28 kDa

研究背景 (Background)

Apoptotic adaptor molecule that recruits caspase-8 or caspase-10 to the activated Fas (CD95) or TNFR-1 receptors. The


产品名称: FADD Rabbit Monoclonal Antibody
产品货号: AMRe85554




resulting aggregate called the death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation. Active caspase-8 initiates the subsequent cascade of caspases mediating apoptosis . Involved in interferon-mediated antiviral immune response, playing a role in the positive regulation of interferon signaling .

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of FADD in Raw264.7 lysates using FADD antibody.



Western blot analysis of FADD in SH-SY5Y lysates using FADD antibody



Immunohistochemistry analysis of paraffin-embedded mouse liver using FADD antibody.High-pressure and temperature

产品名称: FADD Rabbit Monoclonal Antibody
产品货号: AMRe85554



Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .