

产品名称: FADD Rabbit Monoclonal Antibody
产品货号: AMRe85555



产品概述 (Summary)

产品名称 (Production Name)	FADD Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FADD FADD; MORT1; GIG3; Protein FADD; FAS-associated death domain protein;
别名 (Alternative Names)	FAS-associating death domain-containing protein; Growth-inhibiting gene 3 protein; Mediator of receptor induced toxicity
基因 ID (Gene ID)	8772.0
蛋白 ID (SwissProt ID)	Q13158.A synthetic peptide of human FADD

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 28 kDa

产品名称: FADD Rabbit Monoclonal Antibody
产品货号: AMRe85555




研究背景 (Background)

Apoptotic adaptor molecule that recruits caspase-8 or caspase-10 to the activated Fas (CD95) or TNFR-1 receptors. The resulting aggregate called the death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation. Active caspase-8 initiates the subsequent cascade of caspases mediating apoptosis.

研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Western blot analysis of FADD in rat Brain, 3T3 lysates using FADD antibody.

注意事项 (Note)

For research use only .