

产品名称: FAM50A Rabbit Monoclonal Antibody
产品货号: AMRe85557



产品概述 (Summary)

产品名称 (Production Name)	FAM50A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FAM50A
别名 (Alternative Names)	family with sequence similarity 50 member A; 9F; XAP5; HXC26; HXC-26; DXS9928E
基因 ID (Gene ID)	9130.0
蛋白 ID (SwissProt ID)	Q14320.A synthetic peptide of human FAM50A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa; Observed MW: 40 kDa

研究背景 (Background)


产品名称: FAM50A Rabbit Monoclonal Antibody
产品货号: AMRe85557




May be a DNA-binding protein or transcriptional factor.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of FAM50A/XAP in HeLa, CHO-K1, C6, Jurkat lysates using FAM50A antibody.




Western blot analysis of FAM50A/XAP in Smooth Muscle, HCT116, HeLa lysates using FAM50A/XAP antibody.



Immunocytochemistry analysis of FAM50A (green) in HeLa using FAM50A antibody, and DAPI (blue).

产品名称: FAM50A Rabbit Monoclonal Antibody
产品货号: AMRe85557



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using FAM50A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .