

产品名称: FANCB Rabbit Monoclonal Antibody
产品货号: AMRe85559



产品概述 (Summary)

产品名称 (Production Name)	FANCB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FANCB
别名 (Alternative Names)	Fanconi anemia complementation group B; FA2; FAB; FACB; FAAP90; FAAP95
基因 ID (Gene ID)	2187.0
蛋白 ID (SwissProt ID)	Q8NB91.A synthetic peptide of human FANCB

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 98 kDa; Observed MW: 98 kDa

研究背景 (Background)


DNA repair protein required for FANCD2 ubiquitination.

产品名称: FANCB Rabbit Monoclonal Antibody
产品货号: AMRe85559



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of FANCB in 3T3, HeLa lysates using FANCB antibody.

注意事项 (Note)

For research use only .