

产品名称: FES Rabbit Monoclonal Antibody
产品货号: AMRe85563



产品概述 (Summary)

产品名称 (Production Name)	FES Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FES
别名 (Alternative Names)	FPS
基因 ID (Gene ID)	2242.0
蛋白 ID (SwissProt ID)	P07332.A synthetic peptide of human FES

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 93 kDa; Observed MW: 93 kDa

研究背景 (Background)

Tyrosine-protein kinase that acts downstream of cell surface receptors and plays a role in the regulation of the actin


产品名称: FES Rabbit Monoclonal Antibody
产品货号: AMRe85563



cytoskeleton, microtubule assembly, cell attachment and cell spreading. Plays a role in FCER1 (high affinity immunoglobulin epsilon receptor)-mediated signaling in mast cells. Acts down-stream of the activated FCER1 receptor and the mast/stem cell growth factor receptor KIT. Plays a role in the regulation of mast cell degranulation. Plays a role in the regulation of cell differentiation and promotes neurite outgrowth in response to NGF signaling. Plays a role in cell scattering and cell migration in response to HGF-induced activation of EZR. Phosphorylates BCR and down-regulates BCR kinase activity. Phosphorylates HCLS1/HS1, PECAM1, STAT3 and TRIM28.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of FES in Raji lysates using FES antibody.

注意事项 (Note)

For research use only .