

产品名称: FGF1 Rabbit Monoclonal Antibody
产品货号: AMRe85565



产品概述 (Summary)

产品名称 (Production Name)	FGF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGF1 FGF1; FGFA; Fibroblast growth factor 1; FGF-1; Acidic fibroblast growth factor;
别名 (Alternative Names)	aFGF; Endothelial cell growth factor; ECGF; Heparin-binding growth factor 1; HBGF-1
基因 ID (Gene ID)	2246.0
蛋白 ID (SwissProt ID)	P05230.Recombinant protein of human FGF1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 17 kDa; Observed MW: 17 kDa

产品名称: FGF1 Rabbit Monoclonal Antibody
产品货号: AMRe85565




研究背景 (Background)

The heparin-binding fibroblast growth factors play important roles in the regulation of cell survival, cell division, angiogenesis, cell differentiation and cell migration.


研究领域 (Research Area)

TGF-beta signaling pathway, PI3K-Akt signaling pathway, MAPK signaling pathway

图片 (Image Data)



Western blot analysis of FGF1 in mouse heart lysates using FGF1 antibody.



Western blot analysis of FGF1 in C6 lysates using FGF1 antibody.

注意事项 (Note)

For research use only .