

产品名称: FLI1 Rabbit Monoclonal Antibody
产品货号: AMRe85570



产品概述 (Summary)

产品名称 (Production Name)	FLI1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FLI1
别名 (Alternative Names)	FLI1; Friend leukemia integration 1 transcription factor; Proto-oncogene Fli-1; Transcription factor ERGB
基因 ID (Gene ID)	2313.0
蛋白 ID (SwissProt ID)	Q01543.A synthetic peptide of human FLI1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 51 kDa; Observed MW: 51 kDa

研究背景 (Background)


产品名称: FLI1 Rabbit Monoclonal Antibody
产品货号: AMRe85570




FLI1 Sequence-specific transcriptional activator. Recognizes the DNA sequence 5'-C[CA]GGAAGT-3'. Belongs to the ETS family. Can form homodimers or heterodimers with ETV6/TEL1. 2 isoforms of the human protein are produced by alternative splicing.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of FLI1 in Jurkat, C6 lysates using FLI1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FLI1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .