

产品名称: FNTA Rabbit Monoclonal Antibody
产品货号: AMRe85573



产品概述 (Summary)

产品名称 (Production Name)	FNTA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FNTA
别名 (Alternative Names)	FPTA; PTAR2; PGGT1A
基因 ID (Gene ID)	2339.0
蛋白 ID (SwissProt ID)	P49354.A synthetic peptide of human FNTA

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 44 kDa; Observed MW: 44 kDa

研究背景 (Background)

Essential subunit of both the farnesyltransferase and the geranylgeranyltransferase complex. Contributes to the transfer of


产品名称: FNTA Rabbit Monoclonal Antibody
产品货号: AMRe85573



a farnesyl or geranylgeranyl moiety from farnesyl or geranylgeranyl diphosphate to a cysteine at the fourth position from the C-terminus of several proteins having the C-terminal sequence Cys-aliphatic-aliphatic-X. May positively regulate neuromuscular junction development downstream of MUSK via its function in RAC1 prenylation and activation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of FNTA in Jurkat, PC-12, HepG2 lysates using FNTA antibody.

注意事项 (Note)

For research use only .