

产品名称: Folate Binding Protein Rabbit Monoclonal Antibody
产品货号: AMRe85575



产品概述 (Summary)

产品名称 (Production Name)	Folate Binding Protein Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Folate Binding Protein
别名 (Alternative Names)	FOLR1; FOLR; Folate receptor alpha; FR-alpha; Adult folate-binding protein; FBP; Folate receptor 1; Folate receptor; adult; KB cells FBP; Ovarian tumor-associated antigen MOv18
基因 ID (Gene ID)	2348.0
蛋白 ID (SwissProt ID)	P15328.A synthetic peptide of human Folate Binding Protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 30 kDa; Observed MW: 40 kDa

产品名称: Folate Binding Protein Rabbit Monoclonal Antibody
产品货号: AMRe85575




研究背景 (Background)


Binds to folate and reduced folic acid derivatives and mediates delivery of 5-methyltetrahydrofolate to the interior of cells.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Folate Binding Protein in HeLa lysates using Folate Binding Protein antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Folate Binding Protein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .