

产品名称: FRA2 Rabbit Monoclonal Antibody
产品货号: AMRe85581



产品概述 (Summary)

产品名称 (Production Name)	FRA2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FRA2
别名 (Alternative Names)	FOSL2; FRA2; Fos-related antigen 2; FRA-2
基因 ID (Gene ID)	2355.0
蛋白 ID (SwissProt ID)	P15408.A synthetic peptide of human FRA2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 35 kDa; Observed MW: 35-45 kDa

研究背景 (Background)

Controls osteoclast survival and size (By similarity). As a dimer with JUN, activates LIF transcription (By similarity). Activates


产品名称: FRA2 Rabbit Monoclonal Antibody
产品货号: AMRe85581



CEBPB transcription in PGE2-activated osteoblasts.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of FRA2 in K562, HeLa lysates using FRA2 antibody.

注意事项 (Note)

For research use only .