

产品名称: Fragilis Rabbit Monoclonal Antibody
产品货号: AMRe85582



产品概述 (Summary)

产品名称 (Production Name)	Fragilis Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Fragilis
别名 (Alternative Names)	IFITM3; Interferon-induced transmembrane protein 3; Dispanin subfamily A member 2b; DSPA2b; Interferon-inducible protein 1-8U
基因 ID (Gene ID)	10410.0
蛋白 ID (SwissProt ID)	Q01628.A synthetic peptide of human Fragilis

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 15 kDa

研究背景 (Background)


产品名称: Fragilis Rabbit Monoclonal Antibody
产品货号: AMRe85582




IFN-induced antiviral protein that mediates cellular innate immunity to at least three major human pathogens, namely influenza A H1N1 virus, West Nile virus (WNV), and dengue virus (WNV), by inhibiting the early step(s) of replication.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Fragilis in K562, HeLa lysates using Fragilis antibody.



Immunocytochemistry analysis of Fragilis(green) in HeLa using Fragilis antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Fragilis antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Fragilis Rabbit Monoclonal Antibody
产品货号: AMRe85582



注意事项 (Note)

For research use only .