

产品名称: FUBP1 Rabbit Monoclonal Antibody
产品货号: AMRe85586



产品概述 (Summary)

产品名称 (Production Name)	FUBP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FUBP1
别名 (Alternative Names)	DNA helicase V; FBP; FUBP; Fubp1; Fubp4; FUSE-binding protein 1; HDH V
基因 ID (Gene ID)	8880.0
蛋白 ID (SwissProt ID)	Q96AE4.A synthetic peptide of human FUBP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa; Observed MW: 74 kDa

研究背景 (Background)

Regulates MYC expression by binding to a single-stranded far-upstream element upstream of the MYC promoter. May act


产品名称: FUBP1 Rabbit Monoclonal Antibody
产品货号: AMRe85586




both as activator and repressor of transcription.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of FUBP1/FPB in mouse spleen lysates using FUBP1/FPB antibody.




Western blot analysis of FUBP1/FPB in K562, rat Brain, HeLa lysates using FUBP1/FPB antibody



Immunocytochemistry analysis of FUBP1 (green) in 293 using FUBP1 antibody, and DAPI (blue).

产品名称: FUBP1 Rabbit Monoclonal Antibody
产品货号: AMRe85586



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using FUBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .