

**产品名称: FUBP1 Rabbit Monoclonal Antibody**  
**产品货号: AMRe85587**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | FUBP1 Rabbit Monoclonal Antibody       |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,ICC                             |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                            |
|---------------------|--------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                               |
| 修饰 (Modification)   | Unmodified                                                                                 |
| 同种型 (Isotype)       | IgG                                                                                        |
| 克隆 (Clonality)      | Monoclonal                                                                                 |
| 形式 (Form)           | Liquid                                                                                     |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.   |
| 储存溶液 (Buffer)       | Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                      |

## 免疫原信息 (Immunogen)

|                        |                                                                        |
|------------------------|------------------------------------------------------------------------|
| 基因名 (Gene Name)        | FUBP1                                                                  |
| 别名 (Alternative Names) | DNA helicase V; FBP; FUBP; Fubp1; Fubp4; FUSE-binding protein 1; HDH V |
| 基因 ID (Gene ID)        | 8880.0                                                                 |
| 蛋白 ID (SwissProt ID)   | Q96AE4.A synthetic peptide of human FUBP1                              |

## 产品应用 (Application)

|                          |                                               |
|--------------------------|-----------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 68 kDa; Observed MW: 74 kDa    |

## 研究背景 (Background)

Regulates MYC expression by binding to a single-stranded far-upstream element upstream of the MYC promoter. May act


产品名称: FUBP1 Rabbit Monoclonal Antibody  
产品货号: AMRe85587




both as activator and repressor of transcription.

## 研究领域 (Research Area)


## 图片 (Image Data)



Western blot analysis of FUBP1/GBP in mouse brain, mouse spleen lysates using FUBP1/GBP antibody.




Western blot analysis of FUBP1/GBP in Jurkat, C6 lysates using FUBP1/GBP antibody



Immunocytochemistry analysis of FUBP1 (green) in 293 using FUBP1 antibody, and DAPI (blue).

**产品名称: FUBP1 Rabbit Monoclonal Antibody**  
**产品货号: AMRe85587**



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using FUBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

#### **注意事项 (Note)**

For research use only .