

产品名称: GCLC Rabbit Monoclonal Antibody
产品货号: AMRe85601



产品概述 (Summary)

产品名称 (Production Name)	GCLC Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GCLC
别名 (Alternative Names)	GCLC; GLCL; GLCLC; Glutamate--cysteine ligase catalytic subunit; GCS heavy chain; Gamma-ECS; Gamma-glutamylcysteine synthetase
基因 ID (Gene ID)	2729.0
蛋白 ID (SwissProt ID)	P48506.A synthetic peptide of human GCLC

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 73 kDa; Observed MW: 73 kDa

研究背景 (Background)


产品名称: GCLC Rabbit Monoclonal Antibody
产品货号: AMRe85601




Sulfur metabolism; glutathione biosynthesis; glutathione from L-cysteine and L-glutamate: step 1/2.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GCLC in mouse liver lysates using GCLC antibody.



Western blot analysis of GCLC in C6 lysates using GCLC antibody.

注意事项 (Note)

For research use only .