

产品名称: GDF15 Rabbit Monoclonal Antibody
产品货号: AMRe85603



产品概述 (Summary)

产品名称 (Production Name)	GDF15 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GDF15 GDF15; MIC1; PDF; PLAB; PTGFB; Growth/differentiation factor 15; GDF-15;
别名 (Alternative Names)	Macrophage inhibitory cytokine 1; MIC-1; NSAID-activated gene 1 protein; NAG-1; NSAID-regulated gene 1 protein; NRG-1; Placental TGF-beta; Placental bone morphogenetic
基因 ID (Gene ID)	9518.0
蛋白 ID (SwissProt ID)	Q99988.A synthetic peptide of human GDF15

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa; Observed MW: 34 kDa

产品名称: GDF15 Rabbit Monoclonal Antibody
产品货号: AMRe85603




研究背景 (Background)

Highly expressed in placenta, with lower levels in prostate and colon and some expression in kidney.

研究领域 (Research Area)

TGF-beta signaling pathway

图片 (Image Data)



Western blot analysis of GDF15 in HeLa lysates using GDF15 antibody.

注意事项 (Note)

For research use only .