

**产品名称: Gelsolin Rabbit Monoclonal Antibody**  
**产品货号: AMRe85604**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | Gelsolin Rabbit Monoclonal Antibody    |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,ICC                             |
| 种属反应性 (Reactivity)     | Human,Mouse                            |

## 产品性能 (Performance)

|                     |                                                                                            |
|---------------------|--------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                               |
| 修饰 (Modification)   | Unmodified                                                                                 |
| 同种型 (Isotype)       | IgG                                                                                        |
| 克隆 (Clonality)      | Monoclonal                                                                                 |
| 形式 (Form)           | Liquid                                                                                     |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.   |
| 储存溶液 (Buffer)       | Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                      |

## 免疫原信息 (Immunogen)

|                        |                                              |
|------------------------|----------------------------------------------|
| 基因名 (Gene Name)        | Gelsolin                                     |
| 别名 (Alternative Names) | ADF; AGEL; Brevin; Gelsolin; Gsn             |
| 基因 ID (Gene ID)        | 2934.0                                       |
| 蛋白 ID (SwissProt ID)   | P06396.A synthetic peptide of human Gelsolin |

## 产品应用 (Application)

|                          |                                               |
|--------------------------|-----------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 86 kDa; Observed MW: 86 kDa    |

## 研究背景 (Background)

Calcium-regulated, actin-modulating protein that binds to the plus (or barbed) ends of actin monomers or filaments,


**产品名称: Gelsolin Rabbit Monoclonal Antibody**  
**产品货号: AMRe85604**



preventing monomer exchange (end-blocking or capping). It can promote the assembly of monomers into filaments (nucleation) as well as sever filaments already formed. Plays a role in ciliogenesis.

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of Gelsolin in MCF-7, 3T3 lysates using Gelsolin antibody.

## 注意事项 (Note)

For research use only .