

产品名称: GEN1 Rabbit Monoclonal Antibody
产品货号: AMRe85607



产品概述 (Summary)

产品名称 (Production Name)	GEN1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GEN1
别名 (Alternative Names)	Gen
基因 ID (Gene ID)	348654.0
蛋白 ID (SwissProt ID)	Q17RS7.A synthetic peptide of human GEN1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 103 kDa; Observed MW: 103 kDa

研究背景 (Background)

Endonuclease which resolves Holliday junctions (HJs) by the introduction of symmetrically related cuts across the junction


产品名称: GEN1 Rabbit Monoclonal Antibody
产品货号: AMRe85607



point, to produce nicked duplex products in which the nicks can be readily ligated. Four-way DNA intermediates, also known as Holliday junctions, are formed during homologous recombination and DNA repair, and their resolution is necessary for proper chromosome segregation (PubMed:19020614, PubMed:26682650). Cleaves HJs by a nick and counter-nick mechanism involving dual coordinated incisions that lead to the formation of ligatable nicked duplex products. Cleavage of the first strand is rate limiting, while second strand cleavage is rapid. Largely monomeric, dimerizes on the HJ and the first nick occurs upon dimerization at the junction (PubMed:26578604). Efficiently cleaves both single and double HJs contained within large recombination intermediates. Exhibits a weak sequence preference for incision between two G residues that reside in a T-rich region of DNA (PubMed:28049850). Has also endonuclease activity on 5'-flap and replication fork (RF) DNA substrates (PubMed:26578604).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GEN1 in rat Brain, C6, 3T3, HeLa lysates using GEN1 antibody.

注意事项 (Note)

For research use only .