

产品名称: GFAP Rabbit Monoclonal Antibody
产品货号: AMRe85608



产品概述 (Summary)

产品名称 (Production Name)	GFAP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GFAP
别名 (Alternative Names)	GFAP; FLJ45472; cb345; ALXDRD
基因 ID (Gene ID)	2670.0
蛋白 ID (SwissProt ID)	P14136.A synthetic peptide of human GFAP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 50 kDa

研究背景 (Background)

GFAP is commonly used as a marker for intracranial and intraspinal tumors arising from astrocytes. In addition, GFAP


产品名称: GFAP Rabbit Monoclonal Antibody
产品货号: AMRe85608




intermediate filaments are also present in nonmyelin-forming Schwann cells in the peripheral nervous system

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GFAP in rat Brain, Hela lysates using GFAP antibody.



Immunohistochemistry analysis of paraffin-embedded Human Brain using GFAP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .