

产品名称: GSTK1 Rabbit Monoclonal Antibody
产品货号: AMRe85619



产品概述 (Summary)

产品名称 (Production Name)	GSTK1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSTK1
别名 (Alternative Names)	GST; GST class kappa; GST13; Gstk1; hGSTK1
基因 ID (Gene ID)	373156.0
蛋白 ID (SwissProt ID)	Q9Y2Q3.A synthetic peptide of human Glutathione S Transferase kappa 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 25 kDa; Observed MW: 25 kDa

研究背景 (Background)


Significant glutathione conjugating activity is found only with the model substrate, 1-chloro-2,4-dinitrobenzene (CDNB).

产品名称: GSTK1 Rabbit Monoclonal Antibody
产品货号: AMRe85619




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GSTK1 in HeLa lysates using GSTK1 antibody.



Immunocytochemistry analysis of Glutathione S Transferase kappa 1(green) in HeLa using Glutathione S Transferase kappa 1 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .