

**产品名称: Glutathione S Transferase mu 1 Rabbit
Monoclonal Antibody**
产品货号: AMRe85620



产品概述 (Summary)

产品名称 (Production Name)	Glutathione S Transferase mu 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Glutathione S Transferase mu 1
别名 (Alternative Names)	MU; H-B; GST1; GTH4; GTM1; MU-1; GSTM1-1; MGC26563; GSTM1a-1a; GSTM1b-1b
基因 ID (Gene ID)	2944.0
蛋白 ID (SwissProt ID)	P09488.A synthetic peptide of human Glutathione S Transferase mu

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 26 kDa

**产品名称: Glutathione S Transferase mu 1 Rabbit
Monoclonal Antibody**
产品货号: AMRe85620




研究背景 (Background)

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Glutathione S Transferase mu 1 in Jurkat, CHO-K1, HeLa lysates using Glutathione S Transferase mu 1 antibody.

注意事项 (Note)

For research use only .