

产品名称: LGR6 Rabbit Monoclonal Antibody
产品货号: AMRe85628



产品概述 (Summary)

产品名称 (Production Name)	LGR6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LGR6
别名 (Alternative Names)	GPCR; Lgr6
基因 ID (Gene ID)	59352.0
蛋白 ID (SwissProt ID)	Q9HBX8.A synthetic peptide of human GPCR LGR6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 104 kDa; Observed MW: 104 kDa

研究背景 (Background)

Orphan receptor.


产品名称: LGR6 Rabbit Monoclonal Antibody
产品货号: AMRe85628




研究领域 (Research Area)

Wnt signaling pathway


图片 (Image Data)



Western blot analysis of GPCR LGR6 in K562, C6, 3T3, HeLa lysates using GPCR LGR6 antibody.



Immunocytochemistry analysis of LGR6 (green) in HeLa using LGR6 antibody, and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using GPCR LGR6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: LGR6 Rabbit Monoclonal Antibody
产品货号: AMRe85628



注意事项 (Note)

For research use only .