

产品名称: RDC1 Rabbit Monoclonal Antibody
产品货号: AMRe85629



产品概述 (Summary)

产品名称 (Production Name)	RDC1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RDC1
别名 (Alternative Names)	RDC1; CXCR7; RDC-1; CMKOR1; CXC-R7; CXCR-7; GPR159
基因 ID (Gene ID)	57007.0
蛋白 ID (SwissProt ID)	P25106.A synthetic peptide of human GPCR RDC1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 41 kDa; Observed MW: 41 kDa

研究背景 (Background)

Plays a regulatory role in CXCR4-mediated activation of cell surface integrins by CXCL12. Required for heart valve


产品名称: RDC1 Rabbit Monoclonal Antibody
产品货号: AMRe85629



development. Acts as coreceptor with CXCR4 for a restricted number of HIV isolates.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GPCR RDC1/CXCR7 in K562, C6, 3T3, HeLa lysates using RDC1 antibody.

注意事项 (Note)

For research use only .