

产品名称: GPRC5B Rabbit Monoclonal Antibody
产品货号: AMRe85632



产品概述 (Summary)

产品名称 (Production Name)	GPRC5B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPRC5B
别名 (Alternative Names)	RAIG2; RAIG-2
基因 ID (Gene ID)	51704.0
蛋白 ID (SwissProt ID)	Q9NZH0.A synthetic peptide of human GPRC5B

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa; Observed MW: 45 kDa

研究背景 (Background)


G-protein coupled receptor family C group 5 member B

产品名称: GPRC5B Rabbit Monoclonal Antibody
产品货号: AMRe85632




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GPRC5B in A549 lysates using GPRC5B antibody.



Western blot analysis of GPRC5B in rat Brain, 3T3 lysates using GPRC5B antibody.

注意事项 (Note)

For research use only .