

产品名称: GSTA3 Rabbit Monoclonal Antibody
产品货号: AMRe85635



产品概述 (Summary)

产品名称 (Production Name)	GSTA3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSTA3
别名 (Alternative Names)	GTA3; GSTA3-3
基因 ID (Gene ID)	2940.0
蛋白 ID (SwissProt ID)	Q16772.A synthetic peptide of human GSTA3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 25 kDa; Observed MW: 25 kDa

研究背景 (Background)

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Catalyzes


产品名称: GSTA3 Rabbit Monoclonal Antibody
产品货号: AMRe85635



isomerization reactions that contribute to the biosynthesis of steroid hormones. Efficiently catalyze obligatory double-bond isomerizations of delta5-androstene-3,17-dione and delta5-pregnene-3,20-dione, precursors to Tosterone and progesterone, respectively.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GSTA3 in C6, 3T3, HeLa lysates using GSTA3 antibody.

注意事项 (Note)

For research use only .