

产品名称: Histone H1.0 Rabbit Monoclonal Antibody
产品货号: AMRe85648



产品概述 (Summary)

产品名称 (Production Name)	Histone H1.0 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Histone H1.0
别名 (Alternative Names)	Histone H1.0; Histone H1(0); Histone H1.0; N-terminally processed; H1F0; H1FV; Histone H5
基因 ID (Gene ID)	3005.0
蛋白 ID (SwissProt ID)	P07305.A synthetic peptide of human Histone H1.0

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa; Observed MW: 30 kDa

研究背景 (Background)


产品名称: Histone H1.0 Rabbit Monoclonal Antibody
产品货号: AMRe85648




Histone H1.0 is a lysine rich member of the H1 family of linker histones. The H1 family of proteins interacts with linker DNA between nucleosomes and mediates compaction into higher order chromatin. Histones H1 are necessary for the condensation of nucleosome chains into higher-order structures. The H1F0 histones are found in cells that are in terminal stages of differentiation or that have low rates of cell division.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Histone H1.0 in C6, CHO-K1 lysates using Histone H1.0 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Histone H1.0 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .