

产品名称: Histone H4 Rabbit Monoclonal Antibody
产品货号: AMRe85656



产品概述 (Summary)

产品名称 (Production Name)	Histone H4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Histone H4
别名 (Alternative Names)	Histone H4
基因 ID (Gene ID)	121504.0
蛋白 ID (SwissProt ID)	P62805.A synthetic peptide of human Histone H4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 11 kDa; Observed MW: 11 kDa

研究背景 (Background)

The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the


产品名称: Histone H4 Rabbit Monoclonal Antibody
产品货号: AMRe85656




nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a member of the histone H4 family.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Histone H4 in HeLa lysates using Histone H4 antibody.



Immunocytochemistry analysis of Histone H4 (green) in HeLa using Histone H4 antibody, and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using Histone H4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Histone H4 Rabbit Monoclonal Antibody
产品货号: AMRe85656



注意事项 (Note)

For research use only .