

产品名称: HOPX Rabbit Monoclonal Antibody
产品货号: AMRe85669



产品概述 (Summary)

产品名称 (Production Name)	HOPX Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HOPX
别名 (Alternative Names)	Hod; Hop; Ob1; Obl; Hdop; Toto; Cameo; AI848177; AW490897; 1110018K11Rik; 1200015P04Rik; 2300002F06Rik
基因 ID (Gene ID)	74318.0
蛋白 ID (SwissProt ID)	Q8R1H0.Recombinant protein of mouse HOPX

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 8 kDa; Observed MW: 12 kDa

研究背景 (Background)


产品名称: HOPX Rabbit Monoclonal Antibody
产品货号: AMRe85669



Atypical homeodomain protein which does not bind DNA and is required to modulate cardiac growth and development. Acts via its interaction with SRF, thereby modulating the expression of SRF-dependent cardiac-specific genes and cardiac development. Prevents SRF-dependent transcription either by inhibiting SRF binding to DNA or by recruiting histone deacetylase (HDAC) proteins that prevent transcription by SRF. Overexpression causes cardiac hypertrophy (PubMed:12297045, PubMed:12297046). Acts as a co-chaperone for HSPA1A and HSPA1B chaperone proteins and assists in chaperone-mediated protein refolding .

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HOPX/HOD in HeLa lysates using HOPX antibody.

注意事项 (Note)

For research use only .