

产品名称: HOXB4 Rabbit Monoclonal Antibody
产品货号: AMRe85670



产品概述 (Summary)

产品名称 (Production Name)	HOXB4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HOXB4
别名 (Alternative Names)	HOXB4; HOX2F; Homeobox protein Hox-B4; Homeobox protein Hox-2.6; Homeobox protein Hox-2F
基因 ID (Gene ID)	3214.0
蛋白 ID (SwissProt ID)	P17483.A synthetic peptide of human HOXB4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 34 kDa

研究背景 (Background)


产品名称: HOXB4 Rabbit Monoclonal Antibody
产品货号: AMRe85670




Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HOXB4 in K562 lysates using HOXB4 antibody.



Immunocytochemistry analysis of HOXB4 (green) in K562 using HOXB4 antibody, and DAPI (blue).

注意事项 (Note)

For research use only .