

产品名称: HP1 alpha Rabbit Monoclonal Antibody
产品货号: AMRe85671



产品概述 (Summary)

产品名称 (Production Name)	HP1 alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HP1 alpha
别名 (Alternative Names)	CBX5; HP1A; Chromobox protein homolog 5; Antigen p25; Heterochromatin protein 1 homolog alpha; HP1 alpha
基因 ID (Gene ID)	23468.0
蛋白 ID (SwissProt ID)	P45973.A synthetic peptide of human HP1 alpha

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa; Observed MW: 22 kDa

研究背景 (Background)


产品名称: HP1 alpha Rabbit Monoclonal Antibody
产品货号: AMRe85671




Component of heterochromatin that recognizes and binds histone H3 tails methylated at 'Lys-9' (H3K9me), leading to epigenetic repression. In contrast, it is excluded from chromatin when 'Tyr-41' of histone H3 is phosphorylated (H3Y41ph). Can interact with lamin-B receptor (LBR). This interaction can contribute to the association of the heterochromatin with the inner nuclear membrane. Involved in the formation of functional kinetochore through interaction with MIS12 complex proteins.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of HP1 alpha in C6, PC-12 lysates using HP1 alpha antibody.




Western blot analysis of HP1 alpha in 3T3, HeLa lysates using HP1 alpha antibody.

产品名称: HP1 alpha Rabbit Monoclonal Antibody
产品货号: AMRe85671



Immunocytochemistry analysis of HP1 alpha(green) in HeLa using HP1 alpha antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using HP1 alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .