

产品名称: HPRT Rabbit Monoclonal Antibody
产品货号: AMRe85674



产品概述 (Summary)

产品名称 (Production Name)	HPRT Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HPRT
别名 (Alternative Names)	HPRT1; HPRT; Hypoxanthine-guanine phosphoribosyltransferase; HGPRT; HGPRTase
基因 ID (Gene ID)	3251.0
蛋白 ID (SwissProt ID)	P00492.A synthetic peptide of human HPRT

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 25 kDa; Observed MW: 25 kDa

研究背景 (Background)


产品名称: HPRT Rabbit Monoclonal Antibody
产品货号: AMRe85674



Converts guanine to guanosine monophosphate, and hypoxanthine to inosine monophosphate. Transfers the 5-phosphoribosyl group from 5-phosphoribosylpyrophosphate onto the purine. Plays a central role in the generation of purine nucleotides through the purine salvage pathway.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HPRT in Hela, CHO-K1, C6, rat Brain, Jurkat lysates using HPRT antibody.

注意事项 (Note)

For research use only .