

产品名称: Hsp70 Rabbit Monoclonal Antibody
产品货号: AMRe85681



产品概述 (Summary)

产品名称 (Production Name)	Hsp70 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Hsp70
别名 (Alternative Names)	HSP72; HSPA1; HSP70I; HSPA1B; HSP70-1; HSP70-1A.
基因 ID (Gene ID)	3303.0
蛋白 ID (SwissProt ID)	P0DMV8.A synthetic peptide of human Hsp70

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 70 kDa

研究背景 (Background)

In conjunction with other heat shock proteins, this protein stabilizes existing proteins against aggregation and mediates the

产品名称: Hsp70 Rabbit Monoclonal Antibody
产品货号: AMRe85681




folding of newly translated proteins in the cytosol and in organelles. It is also involved in the ubiquitin-proteasome pathway through interaction with the AU-rich element RNA-binding protein 1. The gene is located in the major histocompatibility complex class III region, in a cluster with two closely related genes which encode similar proteins.


研究领域 (Research Area)

MAPK signaling pathway


图片 (Image Data)



Western blot analysis of Hsp70 in K562, rat Brain, C6, 3T3 lysates using Hsp70 antibody.



Immunocytochemistry analysis of Hsp70 (green) in HeLa using Hsp70 antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Hsp70 antibody. High-pressure and temperature

产品名称: Hsp70 Rabbit Monoclonal Antibody
产品货号: AMRe85681



Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .