

产品名称: ICAD Rabbit Monoclonal Antibody
产品货号: AMRe85685



产品概述 (Summary)

产品名称 (Production Name)	ICAD Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ICAD
别名 (Alternative Names)	DFFA; DFF1; DFF45; H13; DNA fragmentation factor subunit alpha; DNA fragmentation factor 45 kDa subunit; DFF-45; Inhibitor of CAD; ICAD
基因 ID (Gene ID)	1676.0
蛋白 ID (SwissProt ID)	O00273.A synthetic peptide of human ICAD

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 37 kDa; Observed MW: 45 kDa

研究背景 (Background)

产品名称: ICAD Rabbit Monoclonal Antibody
产品货号: AMRe85685




Inhibitor of the caspase-activated DNase (DFF40).

研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Western blot analysis of ICAD in rat Brain, MCF-7, C6 lysates using ICAD antibody.

注意事项 (Note)

For research use only .