

产品名称: ICOS Ligand Rabbit Monoclonal Antibody
产品货号: AMRe85687



产品概述 (Summary)

产品名称 (Production Name)	ICOS Ligand Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ICOS Ligand
别名 (Alternative Names)	inducible T-cell co-stimulator ligand; B7H2; GL50; B7-H2; B7RP1; CD275; ICOSL; LICOS; B7RP-1; ICOS-L
基因 ID (Gene ID)	23308.0
蛋白 ID (SwissProt ID)	O75144.Recombinant protein of human ICOS Ligand

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 75 kDa

研究背景 (Background)


产品名称: ICOS Ligand Rabbit Monoclonal Antibody
产品货号: AMRe85687



Ligand for the T-cell-specific cell surface receptor ICOS. Acts as a costimulatory signal for T-cell proliferation and cytokine secretion; induces also B-cell proliferation and differentiation into plasma cells. Could play an important role in mediating local tissue responses to inflammatory conditions, as well as in modulating the secondary immune response by co-stimulating memory T-cell function .

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ICOS Ligand in K562 lysates using ICOS Ligand antibody.

注意事项 (Note)

For research use only .