

产品名称: IDH2 Rabbit Monoclonal Antibody
产品货号: AMRe85688



产品概述 (Summary)

产品名称 (Production Name)	IDH2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IDH2
别名 (Alternative Names)	IDH; IDP; IDHM; IDPM; ICD-M; D2HGA2; mNADP-IDH
基因 ID (Gene ID)	3418.0
蛋白 ID (SwissProt ID)	P48735.A synthetic peptide of human IDH2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 51 kDa; Observed MW: 43 kDa

研究背景 (Background)

Plays a role in intermediary metabolism and energy production. It may tightly associate or interact with the pyruvate


产品名称: IDH2 Rabbit Monoclonal Antibody
产品货号: AMRe85688




dehydrogenase complex.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of IDH2 in rat heart, rat gastric lysates using IDH2 antibody.



Western blot analysis of IDH2 in 3T3 lysates using IDH2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using IDH2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: IDH2 Rabbit Monoclonal Antibody
产品货号: AMRe85688



注意事项 (Note)

For research use only .