

产品名称: IFT88 Rabbit Monoclonal Antibody
产品货号: AMRe85691



产品概述 (Summary)

产品名称 (Production Name)	IFT88 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IFT88
别名 (Alternative Names)	Recessive polycystic kidney disease protein Tg737 homolog; Tetratricopeptide repeat protein 10; TPR repeat protein 10
基因 ID (Gene ID)	8100.0
蛋白 ID (SwissProt ID)	Q13099.A synthetic peptide of human IFT88

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 93 kDa; Observed MW: 93 kDa

研究背景 (Background)


产品名称: IFT88 Rabbit Monoclonal Antibody
产品货号: AMRe85691




Involved in primary cilium biogenesis. Also involved in autophagy since it is required for trafficking of ATG16L and the expansion of the autophagic compartment.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of IFT88 in mouse testis lysates using IFT88 antibody.



Western blot analysis of IFT88 in Jurkat, C6, HeLa lysates using IFT88 antibody.

注意事项 (Note)

For research use only .