

产品名称: IGFBP2 Rabbit Monoclonal Antibody
产品货号: AMRe85694



产品概述 (Summary)

产品名称 (Production Name)	IGFBP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGFBP2
别名 (Alternative Names)	IGFBP2; BP2; IBP2; Insulin-like growth factor-binding protein 2; IBP-2; IGF-binding protein 2; IGFBP-2
基因 ID (Gene ID)	3485.0
蛋白 ID (SwissProt ID)	P18065.Recombinant protein of human IGFBP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 35 kDa; Observed MW: 35 kDa

研究背景 (Background)


产品名称: IGFBP2 Rabbit Monoclonal Antibody
产品货号: AMRe85694



Inhibits IGF-mediated growth and developmental rates. IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of IGFBP2 in T47D lysates using IGFBP2 antibody.

注意事项 (Note)

For research use only .