

产品名称: IKK alpha Rabbit Monoclonal Antibody
产品货号: AMRe85699



产品概述 (Summary)

产品名称 (Production Name)	IKK alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IKK alpha CHUK; IKKA; TCF16; Inhibitor of nuclear factor kappa-B kinase subunit alpha; I-kappa-B kinase alpha; IKK-A; IKK-alpha; IkbKA; IkappaB kinase; Conserved helix-loop-helix ubiquitous kinase; I-kappa-B kinase 1; IKK1; Nuclear factor NF-kappa-B
别名 (Alternative Names)	
基因 ID (Gene ID)	1147.0
蛋白 ID (SwissProt ID)	O15111.A synthetic peptide of human IKK alpha

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 85 kDa; Observed MW: 85 kDa

产品名称: IKK alpha Rabbit Monoclonal Antibody
产品货号: AMRe85699




研究背景 (Background)

Plays an essential role in the NF-kappa-B signaling pathway which is activated by multiple stimuli such as inflammatory cytokines, bacterial or viral products, DNA damages or other cellular stresses. Activation of IKK depends upon phosphorylation at Ser177 and Ser181 in the activation loop of IKK β (Ser176 and Ser180 in IKK α), which causes conformational changes, resulting in kinase activation.


研究领域 (Research Area)

PI3K-Akt signaling pathway, MAPK signaling pathway

图片 (Image Data)



Western blot analysis of IKK alpha in HeLa lysates using IKK alpha antibody.



Western blot analysis of IKK alpha in K562, C6, 3T3, HeLa lysates using IKK alpha antibody.

注意事项 (Note)

For research use only .