

产品名称: Inhibin beta A Rabbit Monoclonal Antibody
产品货号: AMRe85708



产品概述 (Summary)

产品名称 (Production Name)	Inhibin beta A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Inhibin beta A
别名 (Alternative Names)	INHBA; Inhibin beta A chain; Activin beta-A chain; Erythroid differentiation protein; EDF
基因 ID (Gene ID)	3624.0
蛋白 ID (SwissProt ID)	P08476.A synthetic peptide of human Inhibin beta A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 44 kDa

研究背景 (Background)


产品名称: Inhibin beta A Rabbit Monoclonal Antibody
产品货号: AMRe85708



Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary gland. Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival, embryonic axial development or bone growth, depending on their subunit composition.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Inhibin beta A in rat Brain, C6, CHO-K1, HeLa lysates using Inhibin beta A antibody.

注意事项 (Note)

For research use only .