

产品名称: Integrin beta 5 Rabbit Monoclonal Antibody
产品货号: AMRe85712



产品概述 (Summary)

产品名称 (Production Name)	Integrin beta 5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Integrin beta 5
别名 (Alternative Names)	Integrin beta-5
基因 ID (Gene ID)	3693.0
蛋白 ID (SwissProt ID)	P18084.Recombinant protein of human Integrin beta 5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 88 kDa; Observed MW: 88 kDa

研究背景 (Background)

Integrin alpha-V/beta-5 (ITGAV:ITGB5) is a receptor for fibronectin. It recognizes the sequence R-G-D in its ligand.

产品名称: Integrin beta 5 Rabbit Monoclonal Antibody
产品货号: AMRe85712




(Microbial infection) Integrin ITGAV:ITGB5 acts as a receptor for adenovirus type C.

研究领域 (Research Area)

PI3K-Akt signaling pathway

图片 (Image Data)



Western blot analysis of Integrin beta 5 in C6, CHO-K1, HeLa lysates using Integrin beta 5 antibody.

注意事项 (Note)

For research use only .