

产品名称: ISG15 Rabbit Monoclonal Antibody
产品货号: AMRe85717



产品概述 (Summary)

产品名称 (Production Name)	ISG15 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ISG15
别名 (Alternative Names)	IP17; Ubiquitin cross-reactive protein; hUCRP; ISG15; G1P2
基因 ID (Gene ID)	9636.0
蛋白 ID (SwissProt ID)	P05161.A synthetic peptide of human ISG15

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 18 kDa; Observed MW: 18 kDa

研究背景 (Background)

May serve as a trans-acting binding factor directing the association of ligated target proteins to intermediate filaments.


产品名称: ISG15 Rabbit Monoclonal Antibody
产品货号: AMRe85717



May also be involved in autocrine, paracrine and endocrine mechanisms, as in cell-to-cell signaling, possibly partly by inducing IFN-gamma secretion by monocytes and macrophages. Seems to display antiviral activity during viral infections.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ISG15 in C6, 3T3, HeLa lysates using ISG15 antibody.

注意事项 (Note)

For research use only .