

产品名称: Islet 1 Rabbit Monoclonal Antibody
产品货号: AMRe85718



产品概述 (Summary)

产品名称 (Production Name)	Islet 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Islet 1
别名 (Alternative Names)	ISL1; Insulin gene enhancer protein ISL-1; Islet-1
基因 ID (Gene ID)	3670.0
蛋白 ID (SwissProt ID)	P61371.Recombinant protein of human Islet 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 39 kDa; Observed MW: 39 kDa

研究背景 (Background)

ISL1 (ISL1 transcription factor, LIM/homeodomain) is a member of the LIM/homeodomain family of transcription factors. It


产品名称: Islet 1 Rabbit Monoclonal Antibody
产品货号: AMRe85718



binds to the enhancer region of the insulin gene, among others, and may play an important role in regulating insulin gene expression. It is central to the development of pancreatic cell lineages and may also be required for motor neuron generation. Islet-1 expression defines cardiac progenitor cell populations and is required for normal cardiac development and asymmetry.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Islet 1 in K562, rat Brain lysates using Islet 1 antibody.

注意事项 (Note)

For research use only .